

# Evac+Chair® 400H

THE WORLD'S BEST SELLING EVACUATION CHAIR JUST GOT EVEN BETTER!

The Evac+Chair® 400H featuring a weight capacity of 440lbs and single person operation.

**MODEL BENEFITS**

- Reinforced Adjustable Friction Tracks
- 440lbs weight capacity
- 5 Point Safety Harness
- Passenger Armrests and Footrest



**FEATURES**

- 1. Adjustable patient head restraint**  
Limits movement, securing the passenger in place
- 2. Adjustable friction tracks (EVAC+TRAC technology)**  
Adjustable self braking reinforced tracks allow for tension change to control speed of descent to suit most staircases
- 3. Ergonomic adjustable handle**  
For ease of deployment and compact storage
- 4. Passenger armrests**  
For added passenger comfort and safety.
- 5. 4 Point quick release safety harness**  
For increased passenger safety
- 6. Heavy duty wheels**  
Ideal for operation over uneven terrain
- 7. Large lockable rear wheels**  
Can be applied in situations where the chair needs to remain in a stationary position
- 8. Passenger footrest**  
Allows for stable passenger seating position during stairway operation
- 9. Dual position adjustable seat**  
Adjustable seat allows the passenger to choose desired seating position

Dimensions	Height	Width	Depth	Weight	Weight Capacity
Folded	41"	20"	8"	33lbs	440lbs

**BENEFITS**

- Designed in the US
- FDA Registered
- Provides a quick, easy and safe emergency evacuation
- Easy to use
- One person operation
- Lightweight
- No heavy lifting required
- Enables small attendants to easily move passengers three times their size
- Easy to open, ready to use in seconds
- Folds for compact storage
- Exceeds fire safety requirements and emergency response guidelines established to support the Americans with Disabilities Act
- Wall mountable

**INCLUDES**

- **LIFETIME WARRANTY**
- 4 Point Quick Release Harness
- Photoluminescent Signage
- Wall Mounting Brackets
- Dust Cover
- User Guide

**VIDEO LINK:**

[evac-chair.com/watch-video/](http://evac-chair.com/watch-video/)

For further information please contact:

**Evac+Chair Canada**

T 647-946-9277 E sales@evac-chair.com

W www.evac-chair.ca

