



B&D Transfer Seat

The BraunAbility B&D Transfer Seat utilizes the original vehicle driver or front passenger seat to give you up to 180-degree turning power, making transfers in your lowered-floor wheelchair accessible vehicle or full-size van easier. This transfer seat often requires no additional drilling into the vehicle floor and leaves the second-row area open for wheelchair storage and other passenger mobility. The seat can also be driven from any position for optimal comfort.

- **Travel: Horizontal: 17.5" - Vertical: 6.25"**
- **Swivel: 100° (180° swivel option)**
- **Weight Capacity: 500 lbs**
- **Interlock safety device prevents unintentional activation**
- **Manual override system**
- **Controlled with handheld pendant or side-mounted switchbox**
- **Financing available through BraunAbility Finance**





Unoccupied Wheelchair Restraint - UWR

The BraunAbility UWR is the first of its kind – this locking restraint offers a powered release that doesn't add bolts or pins to your manual wheelchair. Simply lock into place in your lowered-floor accessible vehicle or full-size van and transfer into the driver or passenger seats. Your manual chair will remain in an upright position to make your transfer into the chair simpler and safer than ever before.

- Suitable for manual wheelchairs less than 50 lbs
- Rigid and folding wheelchairs
- Folds up compactly when not in use to save cabin space
- Equipped with a manual override lever
- Limited 3-year parts/1-year labor warranty
- Financing available through BraunAbility Finance



braunability.com

800.488.0359

