

Aluminum Rehab Shower Commode Chair

HCPCS:



- Can be positioned over a standard toilet or used as a portable self-contained commode
- White aluminum frame is rust resistant
- Attractive padded seat back and arms
- Seat cut-out provides maximum comfort and convenience
- Height adjustable
- Comes with swing-away footrests

PRODUCTS	Description
NRS185006	Self Propel Shower Commode Chair,1/ea

SPECIFICATIONS	Description
Back Height	19"
Carton Shipping Weight	49 lbs.
Overall Length (with riggings)	42"
Overall Width Open	24.5"
Seat Dimensions	18" (W) x 17" (D)
Seat To Floor Height	19.5"-21.5"
Warranty	Limited Lifetime
Weight	41 lbs.
Weight Capacity	275 lbs

PARTS	Description
NRS185006-02	SeatClipW/ScrewsForNRS185006/7, 1/ea, EA
NRS185006-03	BackInsertForNRS185006/7, 1/ea, EA
NRS185006-01	SeatWithInsertForNRS185006/7, 1/ea, EA
NRS185006-04	ArmpadForNRS185006/7, 1/ea, PR
NRS185006-05	ArmrestBracketLeftForNRS185006/7, 1/ea, EA
NRS185006-06	ArmrestBracketRightForNRS185006/7, 1/ea, EA
NRS185006-08	CasterWithLegForNRS185006/7, 1/ea, EA
NRS185006-09	FrameClipForNRS185006/7, 4/pkg, PKG
NRS185006-10	QuickReleaseAxleForNRS185006Only,1/ea, EA
NRS185006-11	SwingawayFootrestForNRS185006/7, 1/pr, PR
NRS185006-12	FootplateForNRS185006/7, 1/pr, PR
NRS185006-13	AngleAdjFtplateMntW/HdwrNRS185006/7,1/ea, EA
NRS185006-14	DowntubePlasticInsrtsForNRS185006/7,1/ea, EA
NRS185006-15	FootrstHgtAdjKnobW/BoltNRS185006/7,1/ea, EA

PARTS**Description**

NRS185006-16	FootrstPivotSaddleW/BoltNRS185006/7,1/ea, EA
NRS185006-17	FootrestFrameLatchForNRS185006/7, 1/ea, EA
NRS185006-18	WheelLockLeftForNRS185006/7, 1/ea, EA
NRS185006-19	WheelLockLRightForNRS185006/7, 1/ea, EA
NRS185006-20	RearAxleReceiverForNRS185006Only, 1/ea, EA
NRS185006-21	BackLockingKnobForNRS185006/7, 1/pr, PR