

# TRILLIUM TYPE 2 Series

THE TRILLIUM TYPE 2 SERIES IS ONE OF THE MOST COMFORTABLE WALKER ON THE MARKET TODAY.

With the new padded seat cushion and wide back strap pad, users can sit in comfort and walk in style. The Trillium Type 2 comes with an all-new braking system, allowing for effortless handle height adjustment and easy maintenance.

## FEATURES

**CENTRE FOLDING** Folding handle is parallel to folding direction. More ergonomically friendly to fold and lift.

**SEAT** Durable waterproof foam cushion

**BRAKE SYSTEM** Semi-concealed cable. Won't snag on clothes and is easy to maintain

**COLLAPSIBLE BASKET** Foldable & removable

**WHEELS** Lightweight 6" polyurethane tires

**COLORS** Carbon Black, Chinatown Blue, Pink, Shiraz Red, Pearl White, Vivid Green



## OPTIONAL ACCESSORIES

- Cup Holder
- Cane Holder
- Phone Holder (4.7")
- Phone Holder (5.5")
- Slow Down Brakes
- Curb Climber



## SPECIFICATIONS

### MINI

13 lbs\*  
18" Seat Width  
18" Seat Height  
25" Width x 23" Length  
29" – 32" Handle Height  
300 lbs Weight Capacity

### REGULAR

13 lbs\*  
18" Seat Width  
21" Seat Height  
25" Width x 24" Length  
31.5" – 34.5" Handle Height  
300 lbs Weight Capacity

### TALL

14 lbs\*  
18" Seat Width  
24" Seat Height  
25" Width x 24" Length  
34.5" – 37.5" Handle Height  
300 lbs Weight Capacity

\*Net Weight of walker without accessories, backstrap and folding basket