XPRESSO WIDE Series



SPECIFICATIONS

WIDE LOW WIDE

31 lbs* 30 lbs* 20" Seat Width 22" Seat Height 19" Seat Height

24.8" Width x 25.5" Length
32.5" - 37.5" Handle Height
450 lbs Weight Capacity
24.8" Width x 25.5" Length
30.5" - 35.5" Handle Height
450 lbs Weight Capacity

THE XPRESSO WIDE SERIES OFFERS THE MOST STABLE AND DURABLE PLATFORM FOR USERS WITH HEAVY-DUTY REQUIREMENTS.

Made with aerospace-grade aluminum and high strength steel, the Xpresso Wide series is designed with a weight capacity of 450 lbs.

FEATURES

INNOVATIVE DESIGN Roman arch provides support & stability

SEAT Center-folding durable waterproof covered foam cushion

BRAKE SYSTEM Cable-free/ Easy to adjust

COLLAPSIBLE BASKET

Foldable & removable

WHEELS 8in heavy-duty polyurethane tires

FRONT & REAR FORKS

CNC-milled aluminum forks for added strength

NEW PATENTED BACK STRAP Added comfort without compromising visibility

COLORS Red, Silver, Black



OPTIONAL ACCESSORIES

- Cup Holder
- Phone Holder
- Slow Down Brakes
- · Eyeglass Holder
- Cane Holder

*Net Weight of walker without accessories, backstrap and folding basket

